

IRIS 9 CMOS CAMERA

KEY FEATURES

- Small 4.25 µm pixels
- Pair with low magnification
 objectives
- High resolution sensor (9 MP)
- 30 fps imaging
- Simple integration
- Compact form factor
- Programmable scan mode to control camera readout, ideal for light-sheet microscopy

TYPICAL APPLICATIONS

- Light sheet microscopy
- Live cell imaging
- Spatial biology
- Micro-plate readers
- Fluorescence microscopy

RELIABILITY

- Three-year warranty
- Extended warranty available

High-Resolution Imaging CMOS Camera

The Iris family of CMOS cameras are designed to provide high-resolution imaging even at low magnifications, covering fields of view of up to 25 mm. Small pixels across a large array allow for the capture of highly detailed images across a wide area, resulting in a high-throughput solution without sacrificing image quality.

The Iris 9 camera is ideal for light-sheet microscopy, or when needing to image large samples at high resolutions. The Iris 9 offers a 9-megapixel array of small 4.25 μ m pixels, ideal for low-magnification scans of large samples.









Whole Tissue Imaging Dr. Steven Pittler



Light Sheet Microscopy Prof. John Girkin

TELEDYNE

IRIS 9 SPECIFICATIONS

SPECIFICATIONS	Camera Performance	
Sensor	GPixel GSense 5130 scientific CMOS sensor	
Active Array Size	2960 x 2960 (9 megapixel)	
Pixel Area	4.25 μm x 4.25 μm (18.06 μm²)	
Sensor Area	12.61 mm x 12.61 mm (18.06 mm diagonal)	
Peak QE%	>73%	
Spectral Response	200 – 1100 nm	
	Rolling shutter	
Readout Modes	Effective global shutter	
	Programmable scan mode (PCle only)	
Digital Binning	2 x 2	
Linearity	> 99%	
Cooling Options	Air cooled (0 °C @ 30 °C ambient, 0.5 e ⁻ /pixel/second dark current)	
Digital Interfaces	USB 3.0	
Digital interfaces	PCIe	
Lens Interfaces	C-mount	
Mounting Points	4 x 1/4"-20 UNC mounting points	
Camera Weight	0.68 kg, 1.5 lbs	

CAMERA MODES

SPECIFICATIONS	Main Imaging Mode	
Bit Depth	16-bit	
Frame Rate (Full Frame)	30 fps (with PCIe)	
Read Noise	1.5 e ⁻	
Cooling (Air)	0 °C	
Line Time	11.26 µsec/line	
Full Well Capacity	10,000 e-	

TRIGGERING MODE	Function	
INPUT TRIGGER MODES		
Trigger First	Sequence triggered on first rising edge	
Edge Trigger	Each frame in sequence triggered by rising edge	
OUTPUT TRIGGER MODES		
Any Row	Expose signal is high while any row is acquiring data	
First Row	Expose signal is high while first row is acquiring data	
EFFECTIVE GLOBAL SHUTTER TRIGGER MODES		
All Row	Expose out signal high for exposure time, maintains exposure time but drops frame rate	
All Row	Expose out signal high for exposure time - readout time. Keeps frame rate but drops exposure time.	
OUTPUT TRIGGER SIGNALS		
Expose Out (up to four signals), Read Out, Trigger Ready		

TELEDYNE

IRIS 9 QE CURVE



IRIS 9 SPEED TABLE

FRAME RATES (HZ)		
ARRAY SIZE	PCIe INTERFACE	USB 3.0 INTERFACE
2960 x 2960	30	16
2960 x 1500	59	32
2960 x 512	174	94
2960 x 128	695	331

IRIS 9 PROGRAMMABLE SCAN MODE

PROGRAMMABLE SCAN MODE	Function	
SCAN MODES		
Auto	Sequence triggered on first rising edge	
Line Delay	Each frame in sequence triggered by rising edge	
Scan Width	Expose signal is high while any row is acquiring data	
SCAN DIRECTION		
Down	Rolling shutter readout begins at the top of the sensor	
Up	Rolling shutter readout begins at the bottom of the sensor	
Down/Up Alternate	Rolling shutter readout alternates direction after starting at the top of the sensor	

TELEDYNE

IRIS 9 DIMENSIONAL OUTLINES (UNIT: MM)



IRIS 9 ACCESSORIES

ACCESSORIES (INCLUDED)	
PCIe Interface Card	Power supply (12V/10A DC)
PCIe Cable	Quick installation guide
Mini-BNC Trigger Cable	Performance and gain test data

FOR MORE INFORMATION REACH OUT ONLINE:

CONTACT US:	photometrics.com/contact
FOR OEM INQUIRIES:	photometrics.com/oem-page
CONTACT SUPPORT:	photometrics.com/contact/support



Teledyne Photometrics is a registered trademark.

Specifications in this datasheet are subject to change. Refer to the Teledyne Photometrics website for most current specifications. © 2024 Teledyne Photometrics.

Revision Date: 2024 08 20